




Chlorine Bleach Pre-Soak Procedure


WHAT YOU NEED

- Vestec 965
- 20-40L water
- Laundry tub



1 Add 500ml of Vestec 965 to 20-40L of water.

2 Add towels, napkins, linens, etc. Soak for 1-4 hours.
Do not leave in solution longer as it may damage the fabrics.




3 Remove fabrics from solution.

4 Add 500ml of Vestec 995 Antichlor to solution.

5 Use chlorine test papers to check that chlorine is zero-ppm.

6 Solution can now be safely disposed of.

- Vestec 995
- Chlorine Test Paper



Parkside Professional Products Ltd. Niagara Falls, ON, L2H 1J5

1-877-480-8127